

OXYGEN ADMINISTRATION EQUIPMENT

Nasal Cannula

Explain the procedure to the patient

Attach to oxygen source and flush nasal cannula and apply to patient

Non-Rebreather Mask

Explain the procedure to the patient.

Attach to oxygen source and flush tubing and hold one-way valve closed until reservoir is full.

Apply to patient.

Bag-Valve Mask (BVM)

Attach BVM to oxygen source and flush tubing

Extend reservoir and allow to fill with oxygen or allow reservoir bag to fill completely

Apply to patient

Oxygen Flow Rates

Nasal Cannula

1-6 Liters per Minute (LPM)

Non-Rebreather Mask

10-15 Liters per Minute (LPM)

Bag Valve Mask

15 Liters per Minute (LPM)